Fräsgerät AF 30

Bedienungsanleitung (D)

Operation manual (GB)

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Istruzione d'uso (I)

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Made in Switzerland

Nr. 31789

www.nouvag.com

NOUVAG®
Dental and Medical Equipment

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Congratulations on your purchase of a NOUVAG AG product. Thank you for the confidence shown in our products.

Please consult the instruction manual for the use and maintenance of the device in order to ensure that it will function properly and efficiently for many years.

You will find the conformity statement and list of authorized service representatives attached.

- **Before operating, please read instructions carefully!**

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TABLE OF CONTENTS

1	Product description	2
1.1	Use and operation	2
1.2	Technical data	2
1.3	Ambient conditions for operation	2
1.4	Ambient conditions for transport and storage	2
1.5	Device symbols	2
2	Safety instructions	3
3	Supplied equipment	3
4	Device overview	4
5	Commissioning	5
5.1	Device assembly	5
5.2	Milling handpiece and micromotor connection	6
5.3	Modelling table	7
5.4	Drawing pen and measuring tool set	8
5.5	Blockout with 4 blades	9
6	Operation	10
6.1	Switch device on	10
6.2	Setting the rotation speed	10
6.3	Setting the direction of rotation	10
6.4	Using the ON/OFF pedal switch	10
6.5	Without the pedal switch	10
6.6	Switching out the tool tip	11
6.7	Temperature regulation of the blockout	11
7	Maintenance and cleaning	11
8	Faults and fault detection	11
9	Spare part list with order numbers	12
10	Discarding instructions	12

1 Product description

1.1 Use and operation

The AF 30 supports the model analysis, surveying work, construction of attachments, wax scraping and milling, milling of metals in simple milling operations and blocking out, etc.
It has a universal milling unit with milling head and handpiece for milling metals in simple milling operations and a heated wax knife.

1.2 Technical data

Operating voltage: 115V~ / 60 Hz or 230V~ / 50 Hz (see device label)
Power: 130 VA
Fuses: 2 x T2A
Motor speed: 1000 to 30000 rpm
Dimensions (W x H x D): 240 x 335 x 320 mm
Net weight: 10.8 kg
Height adjustment range: 25 mm

1.3 Ambient conditions for operation

Rel. air humidity: max. 80%
Temperature: 10 to 40°C
Air pressure: 800 to 1060 hPa

1.4 Ambient conditions for transport and storage

Rel. air humidity: max. 90%
Temperature: 0 to 60°C
Air pressure: 700 to 1060hPa

1.5 Device symbols

: Caution: refer to enclosed documentation

: EU conformity symbol

: Old electrical and electronic equipment must be disposed separately and may not be included in regular domestic waste.

2 Safety instructions

We are of course very much concerned about your safety and your team's safety. Therefore, it is imperative for you to comply with the following instructions:

- The AF 30 may only be operated by experienced trained personnel!
- Never reach into running or running out drilling or milling machines to prevent risks of injury!
- Never activate the clamping mechanisms as the motor is running or running out! Instruments may get damaged.
- The use of third-party products is the responsibility of the operator!
- Repairs are to be carried out by authorised NOUVAG service technicians only!
- Improper use and repair of the device as well as not complying with the instructions relieve us from any obligations regarding guarantee or other claims!


3 Supplied equipment

Art. No.	Description	Quantity
3028.....	Control unit AF 30	1 unit
1295.....	Micromotor 31	1 unit
1055.....	Milling handpiece.....	1 unit
1888.....	Milling handpiece retaining clip	1 unit
30110.....	Double pivot arm	1 unit
1320.....	Drawing pen	1 unit
1488.....	ON/OFF pedal switch.....	1 unit
1314.....	Modelling table	1 unit
1326.....	Magnetic holding plate	1 unit
1321.....	Blockout Set	1 unit
1343.....	Measuring tool, 3-piece.....	1 unit
30025.....	Stand	1 unit
46503.....	Hexagon wrench SW3	1 unit
23534.....	Hexagon wrench SW4	1 unit
Optional accessories:		
1340.....	Coreformer	1 unit
1315.....	Milling plate, large	1 unit
1316.....	Milling plate, small	1 unit
30087.....	Milling straight edge for magnetic table.....	1 unit
1395.....	Modelling table holder	1 unit
1324.....	Milling tool set, 16-piece.....	1 unit
1092.....	Collet chuck 2.35 mm for drawing pen.....	1 unit

4 Device overview


AF 30 front:

1. Magnetic table
2. Modelling table
3. Magnetic table holder
4. "MAGNETIC TABLE ON / OFF"
5. ON/OFF pedal switch
6. Holder for drawing pen
7. Control knob for motor speed
8. Temperature switch for blackout
9. "HEATING OFF / 1 / 2"
10. Blockout socket "BLOCKOUT"
11. Main switch, device "ON / OFF"
12. Direction switch, motor
13. "FORWARD / REVERSE"
14. Motor socket "MOTOR"
15. Holder for milling tools
16. Stand
17. Halogen lamp
18. Double pivot arm
19. Support for double pivot arm
20. Milling motor
21. Milling handpiece
22. Handpiece retaining clip
23. Drawing pen


AF 30 back:

21. Voltage switch
22. Pedal socket "FOOT SWITCH"
23. Power supply cable "LINE"
24. Motor support

21


22

23

24

5 Commissioning


5.1 Device assembly

		
1) Put AF 30 on a secure base.	2a) With the stand aligned, insert into the stand shaft so that...	2b) ...the flat side of the column faces backwards. Use the hexagon wrench to secure.
		
3) Slide the support for double pivot arm over the stand to the desired working height and use the hexagon wrench to secure (position cam must be aligned upward and backward).	4) Slide the double pivot arm over the stand and snap to support – use butterfly handle to secure.	5) Place the modelling table on the magnetic table.




		
6) Plug the ON/OFF foot switch on the back of the base into the "FOOT SWITCH" socket.	7) Ensure that the operating voltage and mains voltage agree!	8) Connect the mains cable to the mains cable socket "LINE".


5.2 Milling handpiece and micromotor connection

1. Motor connector
2. Micromotor
3. Milling handpiece
4. Actuating sleeve
5. Dowel pin
6. Retaining clip


Assembly:

		
1) Snap the retaining clip to the milling handpiece.	2) Clamp milling handpiece with fixed retaining clip with clamping jaw.	3) Set micromotor in milling handpiece and press until the coupling snaps into place. Use light counter motion to check seat.

		
4) Insert micromotor connector into motor socket "MOTOR".	5) Secure motor cable in clamping slot to prevent accidental damage to the cable.	6) Clamp tool into milling handpiece. Clamp by turning the actuating sleeve. Open collet chuck: Turn actuating sleeve in a clockwise direction until it stops. Close collet chuck: Turn actuating sleeve in a counterclockwise direction until it stops.

5.3 Modelling table

1. Clamping lever
2. Spindle screw
3. Clamping jaws
4. Plate

Clamping the model:

The model can be clamp to the modelling table between three clamping jaws and secured by the spindle screw.

Positioning the modelling table:

The modelling table can be positioned as needed on the magnetic table:

- 1) Switch the magnetic table holder "MAGNETIC TABLE ON / OFF" to the OFF position.
- 2) Move the modelling table to the desired position.
- 3) Switch the magnetic table holder "MAGNETIC TABLE ON / OFF" to the ON position to secure the modelling table.


Adjusting the position of the plate:

The plate can be rotated and pivoted by loosening the clamping lever:

- 1) Loosen clamping cable with sideways motion to the left.
- 2) Turn and pivot the plate to position as desired.
- 3) Draw clamping lever again to the right to secure plate in position.

The built-in positioning unit makes it easy to find the horizontal position.


5.4 Drawing pen and measuring tool set

1. Clamping nut
2. Screw cap
3. Collet key
4. Collet chuck

Measuring tool set, 6-piece:

5. Carbon rod
6. Locating rod 1.5 mm
7. Measuring rod 2.5 mm
8. Measuring rod 3.0 mm
9. Multitransfer rod 3.5 mm
10. Support rod

Clamping the drawing pen:

The drawing pen is used in the milling head clamping device of the double pivot arm. The clamping height can be selected according to the model.

Positioning the drawing pen:

The milling head clamping device and the flexibility of the double pivot arm is used to bring the clamped drawing pen into the desired horizontal and vertical position and to secure it.


Clamping the tool tips:

A 3.00 mm collet chuck is included in the scope of delivery to clamp the accessories. Unscrewing the screw cap allows the collet chuck to come out of the drawing pen receiver. A collet chuck with a 2.35 mm diameter can be ordered for tool tips with 2.35 mm chuck diameter. See spare parts lists with order numbers. Please use the included key to tighten the tool.

5.5 Blockout with 4 blades

The electrically heated wax knife consists of a receiver for the blades and four blades (0°, 2°, 4° and 6°).


1. Receiver
2. Connector for blockout
3. Blade 6°
4. Blade 4°
5. Blade 2°
6. Blade 0°


Assembly:

1) Clamp into receiver.

2) Screw blade with desired cutting angle onto receiver.

3) Clamp drawing pen in double pivot arm.

4) Plug connector into blockout socket "BLOCKOUT".

The device is ready to use!

6 Operation

6.1 Switch device on

The green main switch "ON/OFF" switches the control unit on or off respectively.

→ The switch illuminates when the device is on.

6.2 Setting the rotation speed

The control knob for the motor speed is used to adjust the rotation speed of the motor steplessly from 0 to 30,000 rpm. The number scale indicates the rotation speed in 1000 rpm (e.g. 7.5 corresponds to 7,500 rpm).

If the motor is improperly overstressed for a longer period of time, the built-in overload protection interrupts the power supply and stops the motor. In order to restart the motor, move the main switch to the "OFF" position and wait for at least 60 seconds. The device can then be switched on again and the motor started.

6.3 Setting the direction of rotation

Use the "FORWARD / REVERSE" button to set the motor's direction of rotation.

6.4 Using the ON/OFF pedal switch

The ON/OFF pedal switch can be used to operate the device. This leads to the following operating modes:

<i>Pedal switch not pressed down:</i>	<i>Motor is off.</i>
<i>Pedal switch pressed down:</i>	<i>Motor turns at the speed set by the control knob.</i>

6.5 Without the pedal switch

The device can also be operated without the pedal switch. The "FORWARD / REVERSE" direction switch is used to control the motor:

FORWARD:	<i>Motor turns in a clockwise direction at the speed set by the control knob.</i>
Neutral:	<i>Motor is off.</i>
REVERSE:	<i>Motor turns in a counterclockwise direction at the speed set by the control knob.</i>

Before disconnecting the ON/OFF pedal switch, the "FORWARD / REVERSE" switch must be put in neutral to prevent the motor from turning!

6.6 Switching out the tool tip

- 1) Turn actuating sleeve in a clockwise direction until it stops.
- 2) Remove tool tip from the collet chuck.
- 3) Insert new tool tip.
- 4) Turn actuating sleeve in a counterclockwise direction until it stops.

6.7 Temperature regulation of the blockout

The temperature switch "**HEATING OFF / 1 / 2**" can be used to regulate the temperature:

OFF: Heating element off.

1 / 2: Heating element is heated, where the temperature at level 1 is less than at level 2.

The warm-up time is approx. 1 minute.

7 Maintenance and cleaning

The following important points should be observed for maintaining the device and its accessories:

- The vertical stand must always be lightly greased.
- Do no lubricate the micromotor and handpiece!
- Do not clean the device with compressed air!
- The motor and handpiece must be maintained in accordance with the included description!

Milling dust and swarf must be removed with a brush (included in the scope of delivery) after use. A moist cloth may be used to clean the device and pedal switch.

8 Faults and fault detection

<i>Fault:</i>	<i>Cause:</i>	<i>Elimination:</i>
Device does not function at all	Device not switched on	Move the main switch "ON/OFF" to the "ON" position
	Wrong operating voltage	Check mains voltage
	No connection to the mains	Connect control unit to the mains
Problem with micromotor	Micromotor is overloaded	Switch device off, wait 1 minute, switch device on again
	Contact pins bent	Pull out connector and straighten pins, plug in again
Pedal switch does not function	Pedal switch not connected	Connect pedal switch cable to control unit
Halogen lamp does not illuminate	Light defective	Replace defective light with a new one

If a fault cannot be eliminated, please contact the supplier or an authorised service agent. Addresses can be found on the last page of the operating instructions.

9 Spare part list with order numbers

Accessories	Art. No.
Drawing pen with two collet chucks and keys.....	1320
Collet chuck 2.35 mm for drawing pen	1092
Collet chuck 3.00 mm for drawing pen	1094
Clamping nut for drawing pen	30031
Screw cap	20309
Hexagon wrench SW3	46503
Hexagon wrench SW4	23534
Blockout with 4 blades (0°, 2°, 4° and 6°)	1321
4 blades (0°, 2°, 4° and 6°)	1313
Blade 0°	1354
Blade 2°	1355
Blade 4°	1356
Blade 6°	1357
Measuring tool set	1343
Spare bulb 12 VDC	13230
Retaining clip for milling handpiece	1888

For orders of other parts, please contact our customer service agents who will be glad to assist you.

10 Discarding instructions

When discarding the device, device components and accessories, please comply with the issued statutory regulations.

With regard to the preservation of the environment, old equipment may be returned to the distributor or manufacturer.

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NOUVAG GmbH
NOUVAG USA Inc.